

105.6 Drugs of Abuse in Urine, Multianalyte (powder form)

These SRMs are for verifying the accuracy of methods used to determine cocaine, morphine and codeine and opiate-based substances classed as drugs of abuse in humans. Each SRM consists of multianalyte drugs of abuse. SRM 1511 consists of a mixture of five substances -- morphine, codeine, cocaine metabolite, marijuana metabolite and phencyclidine in freeze-dried urine and is provided as a set of three

[Table 105.5 lists the Single Analyte Drugs of Abuse in Urine.](#)

[Technical Contact: michael.welch@nist.gov](mailto:michael.welch@nist.gov)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

Component in (ng/L)							
SRM	Description	Unit Size	Benzoylecgonine	Morphine	Codeine	THC-9-COOH	Phencyclidine
1511	Multi-Drugs of Abuse	Set of 3 bottles	162	309	288	14.1	23.8
2381	Morphine and Codeine	Set of 4 bottles		134 - 580	130 - 560		